CRP Rapid Test

Test Rápido CRP

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Instructions for Use

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Gebruiksaanweisung

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www.panel.qassay.com

QASSAY CRP Rapid Test Kit

FOR PROFESSIONAL USE
PACKAGING SPECIFICATIONS
10 Tests, 15 Tests, 20 Tests

SF010255, SF010256, SF010257

INTENDED USE

The Qassay CRP Rapid Test kit is intended for the in vitro quantitative determination of CRP in human serum, plasma, and whole blood. It is a fluorescence immunoassay for the diagnosis infection and inflammation and predict future cardiovascular diseases (CVD). For in vitro diagnostic professional use only.

SUMMARY

C-reactive protein (CRP) is synthesized by the liver in response to interleukin-6 and well known as one of the classical acute-phase reactants and as a marker of inflammation. The serum CRP level may rise from a normal level of <5 mg/L to 500 mg/L during the body's general, non-specific response to infectious and other acute inflammatory events; high-sensitivity CRP (hs-CRP) is also emerging as the strongest and most independent predictive risk factor for atherosclerosis and cardiovascular diseases(CVD). For people the diagnosis of inflammatory diseases and CVD assessment cutoffs have been recommended as follows:

Concentrations	Clinical reference
<100 mg/L	Low CVD risk (No Inflammation Situation)
100~300 mg/L	Moderate CVD risk (No Inflammation Situation)
> 3,00 mg/L	High CVD risk (In Inflammation Situation)
> 1000 mg/L	There may be more infections (bacterial infections or viral infections)
10,00~20,00 mg/L	Generally indicates viral infections or mild bacterial infection
20,00~50,00 mg/L	Generally indicates moderate bacterial infection
> 50 mg/L	Generally indicates serious bacterial infection

PRINCIPLE

The Qassay CRP Rapid Test kit is based on fluorescence immunoassay technology. The Qassay CRP Rapid Test Kit uses a sandwich immunodetection method, when sample is added to the sample well of the Test Cartridge, the fluorescence- labeled detector CRP antibody binds to CRP antigen in blood specimen. As the sample mixture migrates on the nitrocellulose matrix of test strip by capillary action, the complexes of detector antibody and CRP are captured to CRP antibody that has been immobilized on test strip. Thus the amount of CRP antigen in blood specimen, the more complexes are accumulated on test strip. Signal intensity of fluorescence of detector antibody reflects amount of CRP captured and Qassay Lateral Flow Reader shows CRP concentration in specimen. The default results unit of Qassay CRP Rapid Test Kit is displayed as XXX mg/L from the reader. The working range and the detection limit of the CRP Test system are 0.50~200.00 mg/L and 0.50 mg/L.

PRECAUTIONS

1. The kits for in vitro diagnostic use only.

2. Do not mix components from different kit lots.

3. Do not use test kit beyond the expiration date.

4. Tests should be applied by professionally trained staff working in certified laboratories at some remove from the patient and clinic at which the sample(s) is taken by qualified medical personnel.

5. The test cartridge should remain in its original sealed pouch until ready to use.

6. Do not use the test cartridge if the pouch is punctured or not well sealed.

7. The Test Cartridge and Meter should be used away from vibration and magnetic field. During normal usage, the Test Cartridge may introduce minute vibration, which should be regarded normal.

8. Use separate clean pipette tips and detector buffer vials for different specimens. The pipette tips and detector buffer vials should be used for one specimen only. Discard after single use.

9. Blood specimens, used test cartridges, pipette tips and detector buffer vials are potentially infectious. Proper laboratory safety techniques, handling and disposal methods should be followed in accordance with standard procedures and relevant regulations observed by microbiological hazard materials.

10. The results should be interpreted by the healthcare professional along with clinical findings and other laboratory test results.

COMPONENTS

Main components (ref. SF010255, SF010256, SF010257):

The product is mainly composed of the following items:

Test cartridge	10 tests	15 tests	20 tests
Detector Buffer	10 units	15 units	20 units
Whole blood sampler	10 units	15 units	20 units
Lancet	10 units	15 units	20 units
Alcohol pad	10 units	15 units	20 units
Biohazard tag	10 units	15 units	20 units
Instructions for use	1 unit	1 unit	1 unit

Note: The components of the different batches of kits are not interchangeable. Materials needed but not provided:

- Qassay Lateral Flow Reader (ref. QIA-DUO-USB)
- Transfer pipette set (10 µL, 100 µL size)
- Sample Collection Containers
- Centrifuge (for Plasma/serum only)

STORAGE AND STABILITY

1. Store the product at 2-30°C in its sealed pouch up to 18 months.

2. Cartridge should be used within 60 minutes after opening the pack.

SAMPLE COLLECTION AND STORAGE

The test can be done with whole blood or serum or plasma.

For whole blood obtained by venipuncture:

- Using standard phlebotomy procedure, collect a venipuncture whole blood specimen using a blood collection tube with suitable anticoagulant (EDTA recommended)

2. It is recommended that specimens should be tested immediately. Do not leave the specimens at room temperature for prolonged periods. If the specimens are not tested immediately, they may be stored at 2°C~8°C.

3. It's not suitable to test the whole blood samples which have been stored at 2°C~8°C for more than 2 days.

For whole blood obtained by fingerpick:

- Wipe the fingertip with the alcohol pad, remove the cap of the lancet and press. Gently squeeze the finger to obtain an adequate amount of blood.
- Gently bring the blood collector close to the drop of blood on your finger to fill it. Make sure that the collector is correctly filled.

For serum or plasma:

- Separate the serum/plasma from blood as soon as possible to avoid hemolysis.
- Test should be performed immediately after the specimens have been collected. Do not leave the specimens at room temperature for prolonged periods. Specimens may be stored at 2°C~8 °C for up to 7 days. For long-term storage, specimens should be kept below -20 °C.

TEST PROCEDURE

Download the app from the Google Play Store or Apple Store and sign in with your credentials. If sign up is needed, please contact your distributor. For specific operation of the instrument, follow the instructions on the reader.

Step 1: Preparation

The test and blood sample should be balanced at room temperature before the test. After the foil bag, take out the test card, and place it on a table. Press on "+" button to start a new test, scan the QR code located in the test and follow the procedure. Select Patient ID, Sample ID, and sample type (if applicable).

Step 2: Sampling

Draw 8.5 µL of whole blood or 5 µL of serum/plasma with a transfer pipette and add it to the buffer tube. If using whole blood sampler open the diluent tube by pulling off the foil at the top and insert the blood collector completely into the diluent liquid.

Step 3: Mixing

Mix well the specimen with buffer for 1 minute by tapping or inverting the tube

Step 4: Loading

Take 75 µL of sample mixture and load it onto the sample well of the Test Cartridge. If using the whole blood sampler, place 2 drops in the sample well.

Step 5: Testing

Press Start Timer button and wait for the 3 minutes reaction time. A timer will be displayed on screen. The Qassay shall be connected to the Mobile App before this timer has expired. Insert the test into the Qassay. Once inserted, the animation on the Qassay App will change and you will be requested to extract the test at a constant speed. After some seconds, the test result will be displayed. If asked to repeat the reading, follow the feedback warning for proper action. If the test is not processed within 1 minute after the countdown concludes, the test will be considered invalid, and a new test must be used for retesting.

QUALITY CONTROL

This product used in conjunction with the Lateral Flow Reader contains internal control for routine quality control requirements. This internal control is carried out every time a patient sample is analyzed. This check indicates whether the test cartridge was inserted and read correctly using the Lateral Flow Reader. In invalid internal control result causes an error message in the analyzer indicating that the test should be repeated.

LIMITATIONS

1. This test has been developed for testing whole blood, human serum and plasma specimen only.

2. The results should be evaluated with all clinical and laboratory data available. If hs-CRP test results do not agree with the clinical evaluation, additional tests should be performed.

3. The false positive results may come from cross-reactions with some similar antibodies in blood; and similar epitopes from non-specific components in blood capturing fluorescence labeled antibodies.

4. EDTA other than anticoagulants (e.g. heparin or citrate) is suggested to use for plasma.

5. Other factors may interfere with the product and may cause erroneous results. These include technical or procedural errors, as well as additional substances in blood specimens.

PRODUCT PERFORMANCE INDEX

Accuracy

A comparison study has been performed using 153 human blood samples, CRP concentrations ranging from 0.54 mg/L to 122.20 mg/L, demonstrated good correlation with a commercially available kit.

Comparison between this test kit and the FincareTM CRP (C-reactive protein) Rapid Quantitative Test was conducted for 153 clinical samples, the Correlation Coefficient is 0.98.

Assay range and Detection Limit

- Assay Range: 0.50~200.00 mg/L
- Detection Limit: 0.50 mg/L

Linearity

A serial concentration of CRP controls of 0.50mg/L, 5.00mg/L, 10.00mg/L, 20.00mg/L, 50.00mg/L, 100.00mg/L were each tested for three times, the Correlation Coefficient (R) is≥0.9900.

Precision

Intra-Run

Within-run precision has been determined by using 10 replicates of four specimens containing 3.00 mg/L and 10.00 mg/L CRP. C.V. is ≤20%.

Inter-Run

Between-run precision has been determined by using 10 replicates for each of three lots using CRP specimen levels at 3.00 mg/L and 10.00 mg/L. C.V. is ≤20%.

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