HbA1C Rapid Test Test Rápido HbA1C Test Rapido della HbA1C HbA1C -Schnelltest

Instructions for Use Instrucciones de Uso Istruzioni per l'Uso Gebrauchsanweisung Οδηγίες χρήσης.

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www.panel.qassay.com

QASSAY HbA1C Rapid Test Kit

FOR PROFESSIONAL USE SFI00255, SFI00256, SFI00257
PACKAGING SPECIFICATIONS
10 Test / Kit, 15 Test / Kit, 20 Test / Kit
INTENDED USE
This product is intended for vitro quantitative determination of Glycated Hemoglobin in human whole blood or peripheral blood samples. This product is used in combination with the Qassay Lateral Flow Reader. Clinically, it is mainly used for auxiliary diagnosis of diabetes and monitoring of blood glucose level.

SUMMARY
Glycosylated hemoglobin is the product of the combination of hemoglobin in red blood cells and blood glucose in human blood. The combination of blood glucose and hemoglobin to produce glycosylated hemoglobin is an irreversible reaction, which is directly proportional to the blood glucose concentration and remains for about 120 days. The glycosylated hemoglobin test can usually reflect the blood glucose control of patients in recent 8-12 weeks, and is not affected by factors such as blood drawing time, fasting and insulin use. Normal reference values: 4.0%-6.0%
Notice: Individual reference range is suggested to be established for each laboratory.

PRINCIPLE
This kit uses immunofluorescence double antibody sandwich method to determine the ratio of glycosylated hemoglobin (HbA1c) in blood. Add the sample into the buffer solution, mix it evenly and drop it into the sample well of the detection card. The glycosylated hemoglobin (HbA1c) and hemoglobin (Hb) in the sample are respectively combined with the fluorescent microsphere labeled HbA1c monoclonal antibody 2 and Hb monoclonal antibody 2 in the fluorescent pad to form an sandwich complex, which diffuse forward along the nitrocellulose membrane under the action of capillary. The complex is captured by HbA1c monoclonal antibody 1 and Hb monoclonal antibody 1 fixed on the nitrocellulose membrane. The concentration of HbA1c antigen and Hb antigen in blood is positively correlated with the fluorescence signal intensity. By the Immunofluorescence Analyzer, the ratio of HbA1c to Hb in the sample is calculated quantitatively according to the ratio of fluorescence signal and the calibration curve of the reagent.

PRECAUTIONS
1. This kit is for in vitro diagnostic use only.
2. Do not mix components from different kits.
3. Do not use test kit beyond the expiration date.
4. Do not use the Test Cartridge if the pouch is punctured or not well sealed.
5. Blood specimens, used Test Cartridges, pipette tips and detector buffer vials should be handled and disposed in accordance with standard procedures and relevant regulations of microbiological hazard materials.
6. The results should be interpreted by the physician along with clinical findings and other laboratory test results.

MATERIALS PROVIDED				
Test cartridge	10 tests	15 tests	20 tests	
Detector Buffer	10 units	15 units	20 units	
Whole blood sampler	10 units	15 units	20 units	
Lancet	10 units	15 units	20 units	
Alcohol pad	10 units	15 units	20 units	
Biohazard bag	10 units	15 units	20 units	
Instructions for use	1 unit	1 unit	1 unit	

Note: The components of the different batches of kits are not interchangeable.

MATERIALS NEEDED BUT NOT PROVIDED
• Qassay Lateral Flow Reader (ref. QIA-DUO-USB)
• Specimen collection containers
• Transfer Pipette Set (10 µL, 100 µL, size)
• Centrifuge (for Plasma/serum only)

STORAGE AND STABILITY
1. Store the detector buffer at 2~30°C. The buffer is stable up to 18 months.
2. Store test kit at 2~30°C, shelf life is up to 18 months.
3. Test Cartridge should be used within 60 minutes after opening the pack.

SAMPLE COLLECTION AND PREPARATION
The test can be performed with whole blood and peripheral blood samples.
For Peripheral Blood:
1. Wipe the fingertip with the alcohol pad, remove the cap of the lancet and press. Gently squeeze the finger to obtain an adequate amount of blood.
2. Gently bring the blood collector close to the drop of blood on your finger to fill it. Make sure that the collector is correctly filled.

For Whole Blood Collected by Venipuncture:
1. Using standard phlebotomy procedure, collect a venipuncture whole blood specimen using a blood collection tube with suitable anticoagulant (EDTA recommended)
2. It is recommended that specimens should be tested immediately. Do not leave the specimens at room temperature for prolonged periods. If the specimens are not tested immediately, they may be stored at 2°C~8°C.
3. It's not suitable to test the whole blood samples storing at 2°C~ 8°C for more than 2 days.
4. Samples with severe lipid blood, hemolysis and microbial contamination cannot be used for the detection of this product. Obviously turbid samples will affect the test results.

TEST PROCEDURE

Download the app from the Google Play Store or Apple Store and sign in with your credentials. If sign up is needed, please contact your distributor. For specific operation of the instrument, follow the instructions on the reader.

Step 1: Preparation
The test and blood sample should be balanced at room temperature before the test. Open the foil bag, take out the test card, and place it on a table. Press "+" button to start a new test, scan the QR code located in the test and follow the procedure. Select Patient ID, Sample ID, and sample type (if applicable).
Step 2: Sampling
Draw 10 µL of whole blood with a transfer pipette and add it to the buffer tube. If using whole blood sample, open the diluent tube by peeling off the foil at the top and insert the blood collector completely into the diluent liquid.
Step 3: Mixing
Mix-well the specimen with buffer 10 times by tapping or inverting the tube and incubate for 1 minute.
Step 4: Loading
Take 75 µL of sample mixture and load it onto the sample well of the Test Cartridge. If using the whole blood sampler, place 2 drops in the sample well.
Step 5: Testing
Press Start Timer button and wait for the 5 minutes reaction time. A timer will be displayed on screen. The Qassay will be connected to the Mobile App before this timer has expired.
Insert the test into the Qassay. Once inserted, the animation on the Qassay App will change and you will be requested to extract the test at a constant speed. After some seconds, the test result will be displayed. If asked to repeat the reading, follow the feedback warning for proper action. If the test is not processed within 1 minute after the countdown concludes, the test will be considered invalid, and a new test must be used for retesting.

QUALITY CONTROL
This product used in conjunction with the Qassay Lateral Flow Reader contains internal control for routine quality requirements. This internal control is performed each time when a patient sample is tested, ensuring correct test insertion and proper reading by the device.

LIMITATIONS
1. This kit is only for the detection of human whole blood samples.
2. Triglycride and bilirubin in the sample will interfere with the test results, and the maximum allowable concentrations are 250g/L and 0.1g/L, respectively. Exceeding the interference concentration and seriously polluting the sample may lead to wrong results.
3. Due to technical reasons, operation errors and other sample factors, test results may be wrong..

PERFORMANCE CHARACTERISTICS
Blank detection limit
The blank detection limit ≤ 4.0%
Accuracy
The relative deviation (Bias) is within ±10.0%.
Linear range
The linear range of determination is between 14.0% and 14.0%, and the linear correlation coefficient (R) is ≥ 0.99.
Repeatability
The coefficient of variation within batches 10.0%.
The coefficient of variation between batches ≤ 10.0%.

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Wuhan Life Origin Biotech Jooer Stock Co., Ltd.