QASSAY Vitamin D Rapid Test Kit

FOR PROFESSIONAL USE

SF150255, SF150256, SF150257,

SF150258, SF150259

PACKAGING SPECIFICATIONS

10 Tests, 15 Tests, 20 Tests

INTENDED USE

The Qassay Vitamin D rapid test kit is intended for *in vitro* quantitative determination of 25-OH-VitD in human serum, plasma, whole blood and peripheral blood samples. This test may help understand the metabolic changes of bone.

PRINCIPLE

This kit uses a competitive method to determine the concentration of 25-hydroxyvitamin D in a sample. The sample was added to the sample diluent, mixed well, and then added dropwise to the test card sampling well, and the 25-hydroxyvitamin D in the sample and the fluorescent microsphere-labeled 25-hydroxyvitamin D antibody in the fluorescent capillary combined to form an antigen-antibody complex. Fluorescent microspheres labeled 25-hydroxyvitamin D monoclonal antibodies that do not form antigen-antibody complexes diffuse forward along the nitrocellulose membrane under capillary action and are captured by BSA protein conjugated to 25-hydroxyvitamin D pre-coated on the nitrocellulose membrane T line. The fluorescent microspheres in the fluorescent pad labeled chicken IgY were captured by the rabbit anti-chicken IgY fixed on the nitrocellulose membrane quality control line. The 25-hydroxyvitamin D in the sample, the less complexes accumulate on the detection line, and the intensity of the fluorescent antibody signal is inversely correlated with the concentration of 25-hydroxyvitamin D in the sample.

PRECAUTIONS

- This kit is for *in vitro* diagnostic use only.
- Do not mix components from different kit lots.
- Do not use test kit beyond the expiration date.
- This product should be applied by professionally trained staff working in certified laboratories at some remove from the patient and clinical at which the sample(s) is taken by qualified medical personnel.
- The test cartridge should remain in its original sealed pouch until ready to use.
- Do not use the test cartridge if the pouch is punctured or not well sealed. Discard after single use.
- Do not smoke, eat, or drink in areas in which specimens or kit reagents are handled.
- Blood specimens, used test cartridges, blood collectors and detector buffer vials are potentially infectious. Proper laboratory safety techniques, handling and disposal methods should be followed in accordance with standard procedures and relevant regulations observed by microbiological hazard materials.
- The results should be interpreted by the physician along with clinical findings and other laboratory test results.
- The test will be applied on a routine basis and not in emergency situations.

MATERIALS PROVIDED

Main components (ref. SF150255, SF150256, SF150257):

The product is mainly composed of the following items:			
Test cartridge	10 tests	15 tests	20 tests
Detector Buffer	10 units	15 units	20 units
Whole blood sampler	10 units	15 units	20 units
Lancet	10 units	15 units	20 units
Alcohol pad	10 units	15 units	20 units
Biohazard bag	10 units	15 units	20 units
Instructions for use	1 unit	1 unit	1 unit

Note: The components of the different batches of kits are not interchangeable.

MATERIALS NEEDED BUT NOT PROVIDED

- Qassay Lateral Flow Reader (ref. Q1A-DUO-USB)
- Centrifuge (for Plasma/Serum only)
- Transfer Pipette Set (10 µL, 100 µL size)
- Specimen collection containers

STORAGE AND STABILITY

- Store the detector buffer at 2~30°C. The buffer is stable up to 18 months.
- Store test kit at 2~30°C, shelf life is up to 18 months.
- Test Cartridge should be used within 60 minutes after opening the pack.

SAMPLE COLLECTION AND PREPARATION

The test can be performed with serum, plasma, whole blood and peripheral blood samples.

For Peripheral Blood:

- Wipe the fingertip with the alcohol pad, remove the cap of the lancet and press. Gently squeeze the finger to obtain an adequate amount of blood.
- Gently bring the blood capillary close to the drop of blood on your finger to fill it. Make sure that the collector is correctly filled.

up to 7 days. For long-term storage, specimens should be kept below -20°C.

TEST PROCEDURE

Download the app from the Google Play Store or Apple Store and sign in with your credentials. If sign up is needed, please contact your distributor. For specific operation of the instrument, follow the instructions on the reader.
Step 1: Preparation
The test and blood sample should be balanced at room temperature before the test. Open the foil bag, take out the test card, and place it on a table. Press on "+" button to start a new test, scan the QR code located in the test and follow the procedure. Select Patient ID, Sample ID, and sample type (if applicable).
Step 2: Sampling
Draw 20 µL of serum, plasma or whole blood with a transfer pipette and add it to the buffer tube. If using whole blood sampler open the diluent tube by peeling off the foil at the top and insert the blood collector completely into the diluent liquid.
Step 3: Mixing
Mix well the specimen with buffer for 10 seconds by tapping or inverting the tube.
Step 4: Loading
Take 80 µL of sample mixture and load it onto the sample well of the Test Cartridge. If using the whole blood sampler, place 2 drops in the sample well.
Step 5: Testing
Press Start Timer button and wait for the 15 minutes reaction time. A timer will be displayed on screen. The Qassay shall be connected to the Mobile App before this timer has expired.

Insert the test into the Qassay. Once inserted, the animation on the Qassay App will change and you will be requested to extract the test at a constant speed. After some seconds, the test result will be displayed. If asked to repeat the reading, follow the feedback warning for proper action. If the test is not processed within 1 minute after the countdown conclusion, the test will be considered invalid, and a new test must be used for retesting.

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QUALITY CONTROL

This product used in conjunction with the Lateral Flow Reader contains internal control for routine quality control requirements. This internal control is carried out every time a patient sample is analyzed. This check indicates whether the test cartridge was inserted and read correctly using the Lateral Flow Reader. An invalid internal control result causes an error message in the analyzer indicating that the test should be repeated.

LIMITATIONS

- This test has been developed for testing human serum, plasma and whole blood specimen only.
- The results of this product should be evaluated with all clinical and laboratory data available. If 25-OH-VitD test results do not agree with the clinical evaluation, additional tests should be performed.
- The false positive results may come from cross-reactions with some similar antibodies in blood; and similar epitopes from non-specific components in blood captured by fluorescent labeled antibodies.
- The false negative results may from some unknown substance blocking epitope adhering antibodies, unstable or degenerated 25-OH-VitD that cannot be identified due to prolonged time and temperature and storage condition of sample and reagent.
- Other factors may interfere with this product and may cause erroneous results. These include technical or procedural errors, as well as additional substances in blood specimens.

PERFORMANCE CHARACTERISTICS

Blank detection limit

The blank detection limit < 5.0 ng/mL.

Accuracy

The relative deviation (Bias) is within ±15.0%.

Linearity

The linear range of determination is between [5.0 ng/mL and 100 ng/mL], and the linear correlation coefficient (R) is ≥ 0.98.

Repeatability

The coefficient of variation within batches: 15.0%.

The coefficient of variation between batches < 20.0%.

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Wuhan Life Origin Biotech Joint Stock Co., Ltd.
SUNGO Europe B.V.